

C17orf53 Rabbit pAb

Catalog No.: A18195

Basic Information

Observed MW

Refer to figures

Calculated MW

70kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

Predicted to enable single-stranded DNA binding activity. Involved in DNA synthesis involved in DNA repair and interstrand cross-link repair. Located in site of DNA damage.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

78995

Swiss Prot

Q8N3J3

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 440-647 of human C17orf53 (NP_076937.2).

Synonyms

MCM8IP; C17orf53

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.