

AHCYL2 Rabbit pAb

Catalog No.: A18204

Basic Information

Observed MW

Refer to figures

Calculated MW

67kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

The protein encoded by this gene acts as a homotetramer and may be involved in the conversion of S-adenosyl-L-homocysteine to L-homocysteine and adenosine. Several transcript variants encoding different isoforms have been found for this gene.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

23382

Swiss Prot

Q96HN2

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human AHCYL2 (NP_056143.1).

Synonyms

IRBIT2; ADOHCYASE3; AHCYL2

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.