

EID1 Rabbit pAb

Catalog No.: A18217

Basic Information

Observed MW

Refer to figures

Calculated MW

21kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

Enables histone acetyltransferase binding activity and histone acetyltransferase regulator activity. Involved in cell differentiation and negative regulation of transcription, DNA-templated. Acts upstream of or within negative regulation of transcription by RNA polymerase II. Located in cytoplasmic ribonucleoprotein granule and nucleoplasm.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

23741

Swiss Prot

Q9Y6B2

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-187 of human EID1 (NP_055150.1).

Synonyms

CRI1; EID-1; RBP21; PTD014; C15orf3; PNAS-22; IRO45620; EID1

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.