

TBC1D21 Rabbit pAb

Catalog No.: A18226 **1 Publications**

Basic Information

Observed MW

Refer to figures

Calculated MW

39kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human, Mouse, Rat

Background

Predicted to enable GTPase activator activity. Predicted to be involved in activation of GTPase activity and intracellular protein transport. Located in extracellular exosome.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

161514

Swiss Prot

Q8IYX1

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-336 of human TBC1D21 (NP_699187.1).

Synonyms

MgcRabGAP; TBC1D21

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.