

SCGB3A1 Rabbit pAb

Catalog No.: A18228

Basic Information

Observed MW

Refer to figures

Calculated MW

10kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

Predicted to be involved in positive regulation of myoblast fusion. Located in extracellular space.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

92304

Swiss Prot

Q96QR1

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SCGB3A1 (NP_443095.2).

Synonyms

HIN1; HIN-1; LU105; UGRP2; PnSP-2; SCGB3A1

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.