

SCGB3A1 Rabbit pAb

Catalog No.: A18228

Basic Information

Observed MW

Refer to figures

Calculated MW

10kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Predicted to be involved in positive regulation of myoblast fusion. Located in extracellular space.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

92304

Swiss Prot

Q96QR1

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

HIN1; HIN-1; LU105; UGRP2; PnSP-2; SCGB3A1

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.