

RIOK1 Rabbit pAb

Catalog No.: A18234

Basic Information

Observed MW

Refer to figures

Calculated MW

66kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

The protein encoded by this gene competes with pICln for inclusion in the protein arginine methyltransferase 5 complex. This complex targets substrates for dimethylation. The encoded protein is essential for the last steps in the maturation of 40S subunits.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

83732

Swiss Prot

Q9BRS2

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human RIOK1 (NP_113668.2).

Synonyms

RIO1; AD034; RRP10; bA288G3.1; RIOK1

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.