

COL14A1 Rabbit pAb

Catalog No.: A18268

Basic Information

Observed MW

Refer to figures

Calculated MW

194kDa

Category

Primary antibody

Applications

IF/ICC

Cross-Reactivity

Human, Mouse

Background

This gene encodes the alpha chain of type XIV collagen, a member of the FACIT (fibril-associated collagens with interrupted triple helices) collagen family. Type XIV collagen interacts with the fibril surface and is involved in the regulation of fibrillogenesis.

Recommended Dilutions

IF/ICC 1:50 - 1:200

Immunogen Information

Gene ID

7373

Swiss Prot

Q05707

Immunogen

A synthetic peptide of human COL14A1.

Synonyms

UND; COL14A1

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.