

# SLR1 Rabbit pAb

Catalog No.: A18329 **2 Publications**

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

65kDa

### Category

Primary antibody

### Applications

WB

### Cross-Reactivity

Oryza sativa

## Background

Probable transcriptional regulator that acts as a repressor of the gibberellin (GA) signaling pathway. Probably acts by participating in large multiprotein complexes that repress transcription of GA-inducible genes. Upon GA application, it is degraded by the proteasome, allowing the GA signaling pathway. In contrast, its overexpression prevents the GA signaling pathway and induces a dwarf phenotype.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

4333860

### Swiss Prot

### Immunogen

Recombinant protein of Oryza sativa SLR1.

### Synonyms

GAI; SLR1; OsGAI; OsJ\_12286

## Contact

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.