

PSKH1 Rabbit pAb

Catalog No.: A18374

Basic Information

Observed MW

Refer to figures

Calculated MW

48kDa

Category

Primary antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

Predicted to enable protein kinase activity. Predicted to act upstream of or within determination of left/right symmetry; heart development; and protein phosphorylation. Located in cytosol and nuclear speck.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

5681

Swiss Prot

P11801

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

PSKH1

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.