

# GPR65 Rabbit pAb

Catalog No.: A18397

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

197kDa

### Category

Primary antibody

### Applications

WB, ELISA

### Cross-Reactivity

Human

## Background

Enables G protein-coupled receptor activity. Involved in several processes, including actin cytoskeleton reorganization; activation of GTPase activity; and positive regulation of stress fiber assembly. Located in plasma membrane.

## Recommended Dilutions

**WB** 1:500 - 1:2000

**ELISA** Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

## Immunogen Information

### Gene ID

8477

### Swiss Prot

Q5VT25

### Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

### Synonyms

TDAG8; hTDAG8; GPR65

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.