

MAST2 Rabbit pAb

Catalog No.: A18425

Basic Information

Observed MW

Refer to figures

Calculated MW

196kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Enables phosphatase binding activity. Predicted to be involved in several processes, including peptidyl-serine phosphorylation; regulation of interleukin-12 production; and spermatid differentiation. Predicted to be located in cytoplasm and plasma membrane. Predicted to be active in microtubule cytoskeleton.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

23139

Swiss Prot

Q6P0Q8

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

MTSSK; MAST205; MAST2

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.