

ANKS1A Rabbit pAb

Catalog No.: A18435

Basic Information

Observed MW

Refer to figures

Calculated MW

123kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Predicted to enable ephrin receptor binding activity. Predicted to be involved in ephrin receptor signaling pathway; neuron remodeling; and substrate-dependent cell migration. Predicted to act upstream of or within negative regulation of ubiquitin-dependent protein catabolic process and regulation of ephrin receptor signaling pathway. Located in cytosol and nucleoplasm.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

23294

Swiss Prot

Q92625

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

ANKS1; ANKS1A

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.