

EPC2 Rabbit pAb

Catalog No.: A18447 **1 Publications**

Basic Information

Observed MW

Refer to figures

Calculated MW

91kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Predicted to contribute to histone acetyltransferase activity. Predicted to be involved in regulation of transcription by RNA polymerase II. Predicted to be part of Piccolo NuA4 histone acetyltransferase complex.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

26122

Swiss Prot

Q52LR7

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

EPC-LIKE; EPC2

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.