

RTFDC1 Rabbit pAb

Catalog No.: A18461

Basic Information

Observed MW

Refer to figures

Calculated MW

34kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Mouse, Rat

Background

Enables DNA binding activity. Involved in cellular response to hydroxyurea and regulation of DNA stability. Located in replication fork.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

51507

Swiss Prot

Q9BY42

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

CDA05; RTFDC1; HSPC164; C20orf43; SHUJUN-3

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.