

ADCY10 Rabbit pAb

Catalog No.: A18473 **1 Publications**

Basic Information

Observed MW

48kDa

Calculated MW

187kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human, Mouse

Background

The protein encoded by this gene belongs to a distinct class of adenylyl cyclases that is soluble and insensitive to G protein or forskolin regulation. Activity of this protein is regulated by bicarbonate. Variation at this gene has been observed in patients with absorptive hypercalciuria. Alternatively spliced transcript variants encoding different isoforms have been observed. There is a pseudogene of this gene on chromosome 6.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

55811

Swiss Prot

Q96PN6

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human ADCY10 (NP_060887.2).

Synonyms

SAC; HCA2; SACI; Sacy; hsAC; HEL-S-7a; ADCY10

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates using ADCY10 Rabbit pAb (A18473) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 180s.