

ARHGAP10 Rabbit pAb

Catalog No.: A18492

Basic Information

Observed MW

Refer to figures

Calculated MW

89kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Mouse

Background

Predicted to enable GTPase activator activity. Predicted to be involved in cytoskeleton organization and negative regulation of apoptotic process. Predicted to be located in perinuclear region of cytoplasm and plasma membrane. Predicted to be active in cytosol.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

79658

Swiss Prot

A1A456

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

GRAF2; PSGAP; PS-GAP; ARHGAP10

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.