

EIL1 Rabbit pAb

Catalog No.: A18498 **1 Publications**

Basic Information

Observed MW

Refer to figures

Calculated MW

69kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Solanum lycopersicum

Background

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

543679

Swiss Prot

Q94FV4

Immunogen

Recombinant protein of Solanum lycopersicum EIL1.

Synonyms

EIN3-like; EIL1

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.