

RPAP2 Rabbit pAb

Catalog No.: A18500

Basic Information

Observed MW

70kDa

Calculated MW

70kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Enables RNA polymerase II CTD heptapeptide repeat phosphatase activity. Involved in dephosphorylation of RNA polymerase II C-terminal domain and snRNA transcription. Located in cytosol and nucleolus. Part of RNA polymerase II, holoenzyme.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

79871

Swiss Prot

Q8IXW5

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Rtr1; C1orf82; RPAP2

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Validation Data



Western blot analysis of lysates from HeLa cells, using RPAP2 Rabbit pAb (A18500) at 1:500 dilution.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 90s.