

GARNL3 Rabbit pAb

Catalog No.: A18508

Basic Information

Observed MW

Refer to figures

Calculated MW

113kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Mouse, Rat

Background

Predicted to enable GTPase activator activity. Predicted to be involved in activation of GTPase activity. Predicted to be active in cytoplasm.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

84253

Swiss Prot

Q5VWW2

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

bA356B19.1; GARNL3

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.