

BTBD10 Rabbit pAb

Catalog No.: A18509

Basic Information

Observed MW

Refer to figures

Calculated MW

54kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Predicted to be involved in negative regulation of neuron death; positive regulation of phosphorylation; and type B pancreatic cell proliferation. Located in fibrillar center and nucleoplasm.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

84280

Swiss Prot

Q9BSF8

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

GMRP1; GMRP-1; BTBD10

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.