

# FBF1 Rabbit pAb

Catalog No.: A18516

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

125kDa

### Category

Primary antibody

### Applications

WB,ELISA

### Cross-Reactivity

Rat

## Background

Involved in apical junction assembly and establishment of epithelial cell polarity. Acts upstream of or within cilium assembly. Located in centriole; centrosome; and spindle pole. Part of ciliary transition fiber. Colocalizes with apical junction complex and keratin filament.

## Recommended Dilutions

**WB** 1:500 - 1:2000

**ELISA** Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

## Immunogen Information

### Gene ID

85302

### Swiss Prot

Q8TES7

### Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

### Synonyms

Alb; FBF-1; FBF1

## Contact

☎ | 400-999-6126

✉ | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

🌐 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.