

ZNF276 Rabbit pAb

Catalog No.: A18523

Basic Information

Observed MW

Refer to figures

Calculated MW

67kDa

Category

Primary antibody

Applications

WB,ELISA

Cross-Reactivity

Mouse, Rat

Background

Enables sequence-specific double-stranded DNA binding activity. Predicted to be involved in regulation of transcription by RNA polymerase II. Located in kinetochore.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

92822

Swiss Prot

Q8N554

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

ZADT; CENPZ; CENP-Z; ZFP276; ZNF477; ZNF276

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.