

ACOT4 Rabbit pAb

Catalog No.: A18535

Basic Information

Observed MW

Refer to figures

Calculated MW

46kDa

Category

Primary antibody

Applications

IHC-P,ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Enables acyl-CoA hydrolase activity and succinyl-CoA hydrolase activity. Involved in carboxylic acid metabolic process; saturated monocarboxylic acid metabolic process; and succinyl-CoA metabolic process. Located in peroxisome.

Recommended Dilutions

IHC-P 1:50 - 1:200

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

122970

Swiss Prot

Q8N9L9

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

PTE1B; PTE2B; PTE-Ib; ACOT4

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.