

TMEM37 Rabbit pAb

Catalog No.: A18545

Basic Information

Observed MW

Refer to figures

Calculated MW

21kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Mouse

Background

Predicted to enable calcium channel activity and voltage-gated ion channel activity. Predicted to be involved in calcium ion transmembrane transport and regulation of ion transmembrane transport. Predicted to be integral component of membrane.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

140738

Swiss Prot

Q8WXS4

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

PR; PR1; TMEM37

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.