

NCKAP5 Rabbit pAb

Catalog No.: A18566

Basic Information

Observed MW

Refer to figures

Calculated MW

209kDa

Category

Primary antibody

Applications

WB,ELISA

Cross-Reactivity

Rat

Background

Predicted to be involved in microtubule bundle formation and microtubule depolymerization.
Predicted to be active in microtubule plus-end.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

344148

Swiss Prot

O14513

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

NAP5; ERIH1; ERIH2; NCKAP5

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.