

HERC1 Rabbit pAb

Catalog No.: A18601

Basic Information

Observed MW

Refer to figures

Calculated MW

532kDa

Category

Primary antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

This gen encodes a member of the HERC protein family. This protein stimulates guanine nucleotide exchange on ARF1 and Rab proteins. This protein may be involved in membrane transport processes.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

8925

Swiss Prot

Q15751

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

p532; p619; MDFPMR; HERC1

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.