

EZR/RDX/MSN Rabbit pAb

Catalog No.: A18647

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

7430/4478/5962

Swiss Prot

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 347-447 of human EZR/RDX/MSNEZR (NP_003370.2).

Synonyms

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.