

LOC101736823 Rabbit pAb

Catalog No.: A18679 **2 Publications**

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Bombyx mori

Background

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

101736823

Swiss Prot

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 251-460 of Bombyx mori LOC101736823.

Synonyms

SQSTM1; LOC101736823

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.