

5-Methylcytosine (5mC) Rabbit pAb

Catalog No.: A18805 **1 Publications**

Basic Information

Observed MW

Calculated MW

Category

Primary antibody

Applications

ELISA,DB

Cross-Reactivity

Species independent

Background

Recommended Dilutions

DB 1:500 - 1:2000

Immunogen Information

Gene ID **Swiss Prot**

Immunogen

5mC

Synonyms

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



The 5mC rabbit polyclonal antibody (A18805) are tested in Dot Blot against 5-methylcytidine(m5C) and 5-methylcytosine(5mC).



Dot-blot analysis of all sorts using 5-Methylcytosine (5mC) Rabbit pAb (A18805)
Line 1 293F cells transfected with TET2 Protein
Line 2 Another 293F cells transfected with TET2 Protein
Line 3 293F cells transfected with GFP-3 Protein
Line 4 DNMT1 knockout (KO) 293F cells