

Nkapl Rabbit pAb

Catalog No.: A18845

Basic Information

Observed MW

Calculated MW

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Mouse

Background

Predicted to enable chromatin binding activity. Predicted to be involved in regulation of gene expression. Predicted to act upstream of or within cell differentiation and spermatogenesis. Predicted to be active in nucleus. Orthologous to human NKAPL (NFKB activating protein like).

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

66707

Swiss Prot

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

4921504I05Rik; Nkapl

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.