

Vrtn Rabbit pAb

Catalog No.: A18847

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Primary antibody

Applications

WB,ELISA

Cross-Reactivity

Mouse

Background

Predicted to be involved in regulation of transcription by RNA polymerase II. Predicted to be active in chromatin. Orthologous to human VRTN (vertebrae development associated).

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

432677

Swiss Prot

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

7420416P09Rik; Vrtn

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.