

N6,2'-O-dimethyladenosine/m6Am Rabbit pAb

Catalog No.: A18871

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Primary antibody

Applications

DB, ELISA

Cross-Reactivity

Human, Mouse, Rat, Other (Wide Range Predicted)

Background

Recommended Dilutions

DB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

Swiss Prot

Immunogen

Chemical compounds corresponding to N6,2'-O-dimethyladenosine/m6Am.

Synonyms

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.