

PTMA Rabbit pAb

Catalog No.: A1956

Basic Information

Observed MW

Refer to figures

Calculated MW

12kDa

Category

Primary antibody

Applications

ELISA

Cross-Reactivity

Human

Background

Enables DNA-binding transcription factor binding activity. Involved in negative regulation of apoptotic process. Located in cytosol and nucleoplasm.

Recommended Dilutions

ELISA

Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

5757

Swiss Prot

P06454

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

TMSA; PTMA

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.