

Fam171a2 Rabbit pAb

Catalog No.: A20246

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Primary antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

Predicted to be located in membrane. Predicted to be integral component of membrane. Orthologous to human FAM171A2 (family with sequence similarity 171 member A2).

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

217219

Swiss Prot

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.