

ANKRD27 Rabbit pAb

Catalog No.: A20263

Basic Information

Observed MW

Refer to figures

Calculated MW

117kDa

Category

Primary antibody

Applications

WB,IHC-P,ELISA

Cross-Reactivity

Human, Rat

Background

Enables guanyl-nucleotide exchange factor activity and small GTPase binding activity. Involved in endocytic recycling and negative regulation of SNARE complex assembly. Acts upstream of or within early endosome to late endosome transport. Located in endosome; lysosome; and plasma membrane. Implicated in eosinophilic esophagitis.

Recommended Dilutions

WB 1:500 - 1:2000

IHC-P 1:50 - 1:200

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

84079

Swiss Prot

Q96NW4

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

VARP; PP12899; ANKRD27

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.