

fimH Rabbit pAb

Catalog No.: A20293

Basic Information

Observed MW

Refer to figures

Calculated MW

31kDa

Category

Primary antibody

Applications

WB,ELISA

Cross-Reactivity

Escherichia coli

Background

Mutants of FimH have been isolated that auto-aggregate. [More information is available at EcoGene: EG10315]. Type 1, or mannose-sensitive, fimbriae in Escherichia coli mediate binding to receptor structures allowing the bacteria to colonize various host tissues. [More information is available at EcoCyc: EG10315].

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

948847

Swiss Prot

P08191

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

ECK4311; pilE; fimH

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.