

MUC5AC Rabbit pAb

Catalog No.: A20312

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Primary antibody

Applications

WB,ELISA

Cross-Reactivity

Mouse

Background

Acts upstream of or within maintenance of lens transparency. Located in extracellular space. Is expressed in nasal cavity epithelium; respiratory system; stomach; and testis. Used to study dry eye syndrome. Human ortholog(s) of this gene implicated in dry eye syndrome. Orthologous to human MUC5AC (mucin 5AC, oligomeric mucus/gel-forming).

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

17833

Swiss Prot

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

MGM; 2210005L13Rik; MUC5AC

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.