

CCDC85C Rabbit pAb

Catalog No.: A20358

Basic Information

Observed MW

Refer to figures

Calculated MW

45kDa

Category

Primary antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

Predicted to be involved in cerebral cortex development. Located in adherens junction.

Recommended Dilutions

WB	1:500 - 1:2000
Elisa	1:500 - 1:2000

Immunogen Information

Gene ID	Swiss Prot
317762	A6NKD9

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

Contact

	400-999-6126
	cn.market@abclonal.com.cn
	www.abclonal.com.cn

Product Information

Source	Isotype	Purification
Rabbit	IgG	Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% sodium azide,40% glycerol,pH7.4.