

# SHROOM1 Rabbit pAb

Catalog No.: A20359

## Basic Information

### Observed MW

100-130kDa

### Calculated MW

91kDa

### Category

Primary antibody

### Applications

WB,IP,ELISA

### Cross-Reactivity

Human

## Background

SHROOM family members play diverse roles in the development of the nervous system and other tissues (Hagens et al., 2006 [PubMed 16615870]).

## Recommended Dilutions

<b>WB</b>	1:500 - 1:3000
<b>IP</b>	1:500 - 1:1000
<b>Elisa</b>	1:500 - 1:2000

## Immunogen Information

### Gene ID

134549

### Swiss Prot

### Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

### Synonyms

APXL2; SHROOM1

## Contact

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.