

RASSF7 Rabbit pAb

Catalog No.: A20361

Basic Information

Observed MW

40kDa

Calculated MW

40kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Rat

Background

Predicted to be involved in apoptotic process; regulation of microtubule cytoskeleton organization; and signal transduction. Located in centriolar satellite.

Recommended Dilutions

WB	1:500 - 1:1000
Elisa	1:500 - 1:2000

Immunogen Information

Gene ID

8045

Swiss Prot

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

HRC1; FAP88; HRAS1; CFAP88; C11orf13; RASSF7

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.