

WDR72 Rabbit pAb

Catalog No.: A20392 **1 Publications**

Basic Information

Observed MW

Refer to figures

Calculated MW

123kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human, Mouse, Rat

Background

This gene encodes a protein with eight WD-40 repeats. Mutations in this gene have been associated with amelogenesis imperfecta hypomaturation type 2A3. Alternative splicing results in multiple transcript variants.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

256764

Swiss Prot

Q3MJ13

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 600-800 of human WDR72 (NP_877435.3).

Synonyms

A12A3; WDR72

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.