

USP18 Rabbit pAb

Catalog No.: A20406

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Primary antibody

Applications

WB,ELISA

Cross-Reactivity

Mouse

Background

Enables ISG15-specific protease activity. Acts upstream of or within response to bacterium and response to stilbenoid. Predicted to be located in intracellular membrane-bounded organelle. Predicted to be active in cytosol and nucleus. Is expressed in axial skeleton; lung; and mandible. Orthologous to several human genes including USP18 (ubiquitin specific peptidase 18).

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

24110

Swiss Prot

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

UBP43; Ubp15; 1110058H21Rik; USP18

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.