

NREP Rabbit pAb

Catalog No.: A20437

Basic Information

Observed MW

Refer to figures

Calculated MW

8kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Predicted to be involved in axon regeneration; regulation of neuron differentiation; and regulation of transforming growth factor beta receptor signaling pathway. Predicted to be active in cytoplasm and nucleus.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

9315

Swiss Prot

Q16612

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

P311; PTZ17; SEZ17; D4S114; C5orf13; PRO1873; NREP

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.