

# $\alpha$ -zein 19KD Rabbit pAb

Catalog No.: A20446 **1 Publications**

## Basic Information

---

### Observed MW

Refer to figures

### Calculated MW

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Zea mays

## Background

---

## Recommended Dilutions

---

WB 1:500 - 1:2000

## Immunogen Information

---

### Gene ID

Swiss Prot

### Immunogen

Recombinant protein of Zea mays  $\alpha$ -zein 19KD.

### Synonyms

## Contact

---

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.