

MT3 Rabbit pAb

Catalog No.: A20541

Basic Information

Observed MW

Refer to figures

Calculated MW

7kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

This gene is a member of the metallothionein family of genes. Proteins encoded by this gene family are low in molecular weight, are cysteine-rich, lack aromatic residues, and bind divalent heavy metal ions. This gene family member displays tissue-specific expression, and contains a threonine insert near its N-terminus and a glutamate-rich hexapeptide insert near its C-terminus relative to the proteins encoded by other gene family members. It plays an important role in zinc and copper homeostasis, and is induced under hypoxic conditions. The encoded protein is a growth inhibitory factor, and reduced levels of the protein are observed in the brains of individuals with some metal-linked neurodegenerative disorders such as Alzheimer's disease.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

4504

Swiss Prot

P25713

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

GIF; GIFB; GRIF; ZnMT3; MT3

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.