

RFWD3 Rabbit pAb

Catalog No.: A20600

Basic Information

Observed MW

Refer to figures

Calculated MW

85kDa

Category

Primary antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

Enables MDM2/MDM4 family protein binding activity; p53 binding activity; and ubiquitin protein ligase activity. Involved in several processes, including DNA metabolic process; regulation of cell cycle phase transition; and response to ionizing radiation. Located in nucleoplasm and site of DNA damage. Colocalizes with site of double-strand break. Implicated in Fanconi anemia.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

55159

Swiss Prot

Q6PCD5

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

FANCW; RNF201; RFWD3

Contact

 | 400-999-6126

 | cn.market@abclonal.com.cn

 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.