

CCNO Rabbit pAb

Catalog No.: A20609

Basic Information

Observed MW

Refer to figures

Calculated MW

38kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

This gene encodes a member of the cyclin protein family, and the encoded protein is involved in regulation of the cell cycle. Disruption of this gene is associated with primary ciliary dyskinesia-19. Alternative splicing results in multiple transcript variants. This gene, which has a previous symbol of UNG2, was erroneously identified as a uracil DNA glycosylase in PubMed ID: 2001396. A later publication, PubMed ID: 8419333, identified this gene's product as a cyclin protein family member. The UNG2 symbol is also used as a specific protein isoform name for the UNG gene (GeneID 7374), so confusion exists in the scientific literature and in some databases for these two genes.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

10309

Swiss Prot

P22674

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

CCNU; UDG2; CILD29; CCNO

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.