

# Xa21 Rabbit pAb

Catalog No.: A20643

## Basic Information

### Observed MW

108kDa

### Calculated MW

111kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Oryza sativa

## Background

Bacterial Blight resistance 21

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

9271075

### Swiss Prot

Q40640

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 898-998 of Oryza sativa Xa21. (Q2R2K5).

### Synonyms

Xa21

## Contact

 | 400-999-6126

 | [cn.market@abclonal.com.cn](mailto:cn.market@abclonal.com.cn)

 | [www.abclonal.com.cn](http://www.abclonal.com.cn)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.

## Validation Data

---



Western blot analysis of extracts of various tissues from the japonica rice (*Oryza sativa* L.) variety Zhonghua 11, using Xa21 Rabbit pAb (A20643) at 1:1000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Enhanced Kit (RM00021).

Exposure time: 60s.