

PDZK1IP1 Rabbit pAb

Catalog No.: A20645

Basic Information

Observed MW

Refer to figures

Calculated MW

12kDa

Category

Primary antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Located in extracellular exosome.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

10158

Swiss Prot

Q13113

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

DD96; SPAP; MAP17; PDZK1IP1

Contact

☎ | 400-999-6126

✉ | cn.market@abclonal.com.cn

🌐 | www.abclonal.com.cn

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.